



Standard Instrument Set for Dynamic and Condylar Screws and Plates in Graphic Case

Listing: T-233 Standard Instrument Set for Dynamic Hip and Condylar Screws and Plates



T-234 Tray for Instrument Set, Dynamic Hip and Condylar Screws and Plates (Steel)



19E-100.01 Standard Container, silver, 580 x 280 x 100 mm (recommended)



Standard Instrument Set for Dynamic and Condylar Screws and Plates in Graphic Case

- Listing:** **T-233** **Standard Instrument Set for Dynamic Hip and Condylar Screws and Plates**
- T-234 Tray for Instrument Set, Dynamic Hip and Condylar Screws and Plates (Steel)

Cat. No.:		Pcs.:
2-300	Calibrated Threaded Guide Pin Dia. 2.5mm x 230 mm Length, Threaded Length 10 mm	10
2-302	Drill Guide adjustable to 125° / 130° / 135° / 140° / 145° / 150°	1
2-304	Measuring Gauge	1
2-306	Triple Reamer complete for Dynamic Hip Plate System	1
2-308	Tap	1
2-310	Guide Sleeve, short, for tapping	1
2-312	Guide Sleeve, long, for screwing in of Lag Screws	1
2-314	Guide Shaft for Insertion of Coupling Screws	1
2-316	Coupling Screw for Femoral Dynamic Hip Screws	1
2-317	Guide Shaft	1
2-318	Wrench	1
2-320	Impactor	1
2-322	Handle	1
2-486	Drill Guide 95°	1
2-488	Triple Reamer complete for Dynamic Condylar System	1

Dynamic Hip Screw Standard Implant Set in Graphic Case

Listing: T-2400 **Dynamic Hip Screw Standard Implant Set**
(Steel Plates in Mat. DIN ISO 5832-1)



T-2410 Tray for Dynamic Hip Screw Standard Implant Set (Steel)



19E-100.01 Standard Container, silver, 580 x 280 x 100 mm
(recommended)

Dynamic Hip Screw Standard Implant Set in Graphic Case



**Listing: T-2400 Dynamic Hip Screw Standard Implant Set
(Steel Plates in Mat. DIN ISO 5832-1)**

T-2410 Tray for Dynamic Hip Screw Standard Implant Set (Steel)

Cat. No.:	Dynamic Lag Screw, Thread Length 22 mm	Pcs.:
2-406	Length 65 mm, Steel	1
2-408	Length 70 mm, Steel	1
2-410	Length 75 mm, Steel	1
2-412	Length 80 mm, Steel	1
2-414	Length 85 mm, Steel	1
2-416	Length 90 mm, Steel	1
2-418	Length 95 mm, Steel	1
2-420	Length 100 mm, Steel	1
2-422	Length 105 mm, Steel	1
2-424	Length 110 mm, Steel	1
2-426	Length 115 mm, Steel	1
<hr/>		
2-440	Compression Screw	3
<hr/>		
Cat. No.:	Dynamic Hip Screw Plates 135°	
2-352	4 Holes,Length 78 mm, Steel	3
2-354	5 Holes,Length 94 mm, Steel	2
2-356	6 Holes,Length 110 mm, Steel	2
<hr/>		
Cat. No.:	Dynamic Hip Screw Plates 150°	
2-382	4 Holes,Length 78 mm, Steel	1
2-384	5 Holes,Length 94 mm, Steel	1
2-386	6 Holes,Length 110 mm, Steel	1

Dynamic Condylar Standard Implant Set in Graphic Case

Listing: T-2500 Dynamic Condylar Standard Implant Set
(Steel Plates in Mat. DIN ISO 5832-1)



T-2510 Tray for Dynamic Condylar Standard Implant Set (Steel)



19E-100.01 Standard Container, silver, 580 x 280 x 100 mm
(recommended)


Dynamic Condylar Standard Implant Set in Graphic Case

- Listing:** **T-2500** **Dynamic Condylar Standard Implant Set
(Steel Plates in Mat. DIN ISO 5832-1)**
- T-2510 Tray for Dynamic Condylar Standard Implant Set (Steel)

Cat. No.:	Dynamic Lag Screw, Thread Length 22 mm	Pcs.:
2-400	Length 50 mm, Steel	2
2-402	Length 55 mm, Steel	2
2-404	Length 60 mm, Steel	2
2-406	Length 65 mm, Steel	2
2-408	Length 70 mm, Steel	2
2-410	Length 75 mm, Steel	1
<hr/>		
2-440	Compression Screw	3
<hr/>		
Cat. No.:	Dynamic Condylar Screw Plates 95°	
2-450	6 Holes, Length 100 mm, Steel	2
2-452	8 Holes, Length 130 mm, Steel	2
2-454	10 Holes, Length 163 mm, Steel	2
2-456	12 Holes, Length 198 mm, Steel	1
2-458	14 Holes, Length 225 mm, Steel	1
2-460	16 Holes, Length 260 mm, Steel	1

Broad and Narrow Plates Set in Graphic Case

**Listing: T-330 Broad and Narrow Plates Set
(Steel Plates in Mat. DIN ISO 5832-1)**



T-331 Tray for Broad and Narrow Plates Set



19E-100.01 Standard Container, silver, 580 x 280 x 100 mm
(recommended)

Broad and Narrow Plates Set in Graphic Case

**Listing: T-330 Broad and Narrow Plates Set
(Steel Plates in Mat. DIN ISO 5832-1)**

T-331 Tray for Broad and Narrow Plates Set (Steel)

		Pcs.
Semi-tubular Plates		
3-674	4.5 / 4 Hole, Length 71 mm Steel	2
3-676	4.5, 5 / 5 Hole, Length 87 mm Steel	2
3-678	4.5 / 6 Hole, Length 103 mm Steel	2
3-680	4.5 / 7 Hole, Length 119 mm Steel	1
Narrow Dynamic Compression Plates		
3-604	4.5 / 4 Hole, Length 71 mm Steel	2
3-606	4.5 / 5 Hole, Length 87 mm Steel	2
3-608	4.5 / 6 Hole, Length 103 mm Steel	4
3-610	4.5 / 7 Hole, Length 119 mm Steel	2
3-612	4.5 / 8 Hole, Length 135 mm Steel	2
3-614	4.5 / 9 Hole, Length 151 mm Steel	2
3-616	4.5 / 10 Hole, Length 167 mm Steel	1
3-620	4.5 / 12 Hole, Length 199 mm Steel	1
Broad Dynamic Compression Plates		
3-640	4.5 / 6 Hole, Length 103 mm Steel	2
3-642	4.5 / 7 Hole, Length 119 mm Steel	2
3-644	4.5 / 8 Hole, Length 135 mm Steel	2
3-646	4.5 / 9 Hole, Length 151 mm Steel	1
3-648	4.5 / 10 Hole, Length 167 mm Steel	1
3-652	4.5 / 12 Hole, Length 199 mm Steel	1
3-656	4.5 / 14 Hole, Length 231 mm Steel	1
Spoon Plates		
3-692	4.5 / 5 Hole, Length 100 mm Steel	1
3-694	4.5 / 6 Hole, Length 120 mm Steel	1
T-Plates		
3-702	4.5 / 4 Hole, Length 84 mm Steel	1
3-706	4.5 / 6 Hole, Length 116 mm Steel	1
3-708	4.5 / 8 Hole, Length 148 mm Steel	1
T-Buttress Plate		
3-714	4.5 / 4 Hole, Length 81 mm Steel	1
L-Buttress Plate		
3-720	4.5 / 4 Hole, Length 85 mm left Steel	1
3-722	4.5 / 4 Hole, Length 85 mm right Steel	1

Steel Screw Set in Graphic Case

Listing: T-336 **Screw Set**
(Steel Screws in Mat. DIN ISO 5832-1)



T-337 Tray for Screw Set



19E-191.01 1/2 Standard Container, silver, 285 x 280 x 135 mm
(recommended)

Steel Screw Set in Graphic Case

**Listing: T-336 Screw Plates Set
(Steel Screws in Mat. DIN ISO 5832-1)**

T-337 Tray for Screw Set

Cat. No. Dia 4.5 mm Cortex Screws. Pcs.			Cat. No. Dia 6.5 mm/16 Cortex Screws. Pcs.		
3-200	Length 14 mm	6	3-396	Length 65 mm	2
3-202	Length 16 mm	6	3-398	Length 70 mm	2
3-204	Length 18 mm	6	3-400	Length 75 mm	2
3-206	Length 20 mm	6	3-402	Length 80 mm	2
3-208	Length 22 mm	6	3-404	Length 85 mm	2
3-210	Length 24 mm	6	3-406	Length 90 mm	2
3-212	Length 26 mm	7	3-408	Length 95 mm	2
3-214	Length 28 mm	7	3-410	Length 100 mm	2
3-216	Length 30 mm	14	3-412	Length 105 mm	2
3-218	Length 32 mm	14	3-414	Length 110 mm	2
3-220	Length 34 mm	14			
3-222	Length 36 mm	14	Cat. No. Dia 6.5 mm/16 Cortex Screws. Pcs.		
3-224	Length 38 mm	7	3-436	Length 45 mm	3
3-226	Length 40 mm	14	3-438	Length 50 mm	3
3-228	Length 42 mm	6	3-440	Length 55 mm	3
3-230	Length 44 mm	6	3-442	Length 60 mm	3
3-232	Length 46 mm	6	3-444	Length 65 mm	3
3-234	Length 48 mm	3	3-446	Length 70 mm	3
3-238	Length 52 mm	3	3-448	Length 75 mm	3
3-242	Length 56 mm	3	3-450	Length 80 mm	3
3-246	Length 60 mm	3	3-452	Length 85 mm	3
3-248	Length 65 mm	3	3-454	Length 90 mm	3
3-250	Length 70 mm	3	3-456	Length 95 mm	3
			3-458	Length 100 mm	3
			3-460	Length 105 mm	3
			3-462	Length 110 mm	3
Cat. No. Dia 4.5 mm Malleolar Screws. Pcs.			Cat. No. Dia 6.5 mm/16 mm Cancellous Screws. Pcs.		
3-352	Length 25/12 mm	3	3-480	fully threaded Length 25 mm	4
3-354	Length 30/15 mm	3	3-482	Length 30 mm	4
3-356	Length 35/18 mm	3	3-484	Length 35 mm	4
3-358	Length 40/20 mm	3	3-486	Length 40 mm	4
3-360	Length 45/22 mm	3	3-488	Length 45 mm	2
3-362	Length 50/24 mm	3	3-490	Length 50 mm	2
3-364	Length 55/26 mm	3	3-492	Length 55 mm	2
3-366	Length 60/28 mm	3	3-494	Length 60 mm	2
3-368	Length 65/30 mm	3	3-496	Length 65 mm	2
3-370	Length 70/32 mm	3	3-498	Length 70 mm	2
			3-500	Length 75 mm	2
Cat. No. Dia 6.5 mm/16 mm Cancellous Screws. Pcs.			3-502	Length 80 mm	2
3-382	Length 30 mm	2	3-504	Length 85 mm	2
3-384	Length 35 mm	2	3-506	Length 90 mm	2
3-386	Length 40 mm	2	3-508	Length 95 mm	2
3-388	Length 45 mm	2	3-510	Length 100 mm	2
3-390	Length 50 mm	2	3-512	Length 105 mm	2
3-392	Length 55 mm	2	3-514	Length 110 mm	2
3-394	Length 60 mm	2			
			Cat. No. Dia 6.5 mm/16 mm Cancellous Screws. Pcs.		
			3-522	Washer for Dia. 6.5 mm Cancellous Screws	6
			4-184	Forceps	1

Titanium Screw Set in Graphic Case

Listing: T-3360 **Screw Set**
(Titanium Screws in Mat. DIN ISO 5832-3)



T-337 Tray for Screw Set



19E-191.01 1/2 Standard Container, silver, 285 x 280 x 135 mm
(recommended)

Titanium Screw Set in Graphic Case


Listing: T-3360 Screw Set
(Titanium Screws in Mat. DIN ISO 5832-3)

T-337 Tray for Screw Set

Cat. No.	Dia	Length	Pcs.	Cat. No.	Dia	Length	Pcs.
Cat. No. Dia 4.5 mm Cortex Screws.				Cat. No. Dia 6.5 mm/16 Cortex Screws.			
3-20000	4.5 mm	14 mm	6	3-39600	6.5 mm/16	65 mm	2
3-20200	4.5 mm	16 mm	6	3-39800	6.5 mm/16	70 mm	2
3-20400	4.5 mm	18 mm	6	3-40000	6.5 mm/16	75 mm	2
3-20600	4.5 mm	20 mm	6	3-40200	6.5 mm/16	80 mm	2
3-20800	4.5 mm	22 mm	6	3-40400	6.5 mm/16	85 mm	2
3-21000	4.5 mm	24 mm	6	3-40600	6.5 mm/16	90 mm	2
3-21200	4.5 mm	26 mm	7	3-40800	6.5 mm/16	95 mm	2
3-21400	4.5 mm	28 mm	7	3-41000	6.5 mm/16	100 mm	2
3-21600	4.5 mm	30 mm	14	3-41200	6.5 mm/16	105 mm	2
3-21800	4.5 mm	32 mm	14	3-41400	6.5 mm/16	110 mm	2
3-22000	4.5 mm	34 mm	14	Cat. No. Dia 6.5 mm/16 mm Cancellous Screws.			
3-22200	4.5 mm	36 mm	14	3-43600	6.5 mm/16	45 mm	3
3-22400	4.5 mm	38 mm	7	3-43800	6.5 mm/16	50 mm	3
3-22600	4.5 mm	40 mm	14	3-44000	6.5 mm/16	55 mm	3
3-22800	4.5 mm	42 mm	6	3-44200	6.5 mm/16	60 mm	3
3-23000	4.5 mm	44 mm	6	3-44400	6.5 mm/16	65 mm	3
3-23200	4.5 mm	46 mm	6	3-44600	6.5 mm/16	70 mm	3
3-23400	4.5 mm	48 mm	3	3-44800	6.5 mm/16	75 mm	3
3-23800	4.5 mm	52 mm	3	3-45000	6.5 mm/16	80 mm	3
3-24200	4.5 mm	56 mm	3	3-45200	6.5 mm/16	85 mm	3
3-24600	4.5 mm	60 mm	3	3-45400	6.5 mm/16	90 mm	3
3-24800	4.5 mm	65 mm	3	3-45600	6.5 mm/16	95 mm	3
3-25000	4.5 mm	70 mm	3	3-45800	6.5 mm/16	100 mm	3
Cat. No. Dia 4.5 mm Malleolar Screws.				Cat. No. Dia 6.5 mm/16 mm Cancellous Screws.			
3-35200	4.5 mm	25/12 mm	3	3-48000	6.5 mm/16	25 mm	4
3-35400	4.5 mm	30/15 mm	3	3-48200	6.5 mm/16	30 mm	4
3-35600	4.5 mm	35/18 mm	3	3-48400	6.5 mm/16	35 mm	4
3-35800	4.5 mm	40/20 mm	3	3-48600	6.5 mm/16	40 mm	4
3-36000	4.5 mm	45/22 mm	3	3-48800	6.5 mm/16	45 mm	2
3-36200	4.5 mm	50/24 mm	3	3-49000	6.5 mm/16	50 mm	2
3-36400	4.5 mm	55/26 mm	3	3-49200	6.5 mm/16	55 mm	2
3-36600	4.5 mm	60/28 mm	3	3-49400	6.5 mm/16	60 mm	2
3-36800	4.5 mm	65/30 mm	3	3-49600	6.5 mm/16	65 mm	2
3-37000	4.5 mm	70/32 mm	3	3-49800	6.5 mm/16	70 mm	2
Cat. No. Dia 6.5 mm/16 mm Cancellous Screws.				Cat. No. Dia 6.5 mm/16 mm Cancellous Screws.			
3-38200	6.5 mm/16	30 mm	2	3-50000	6.5 mm/16	75 mm	2
3-38400	6.5 mm/16	35 mm	2	3-50200	6.5 mm/16	80 mm	2
3-38600	6.5 mm/16	40 mm	2	3-50400	6.5 mm/16	85 mm	2
3-38800	6.5 mm/16	45 mm	2	3-50600	6.5 mm/16	90 mm	2
3-39000	6.5 mm/16	50 mm	2	3-50800	6.5 mm/16	95 mm	2
3-39200	6.5 mm/16	55 mm	2	3-51000	6.5 mm/16	100 mm	2
3-39400	6.5 mm/16	60 mm	2	3-51200	6.5 mm/16	105 mm	2
				Cat. No. Dia 6.5 mm/16 mm Cancellous Screws.			
				3-52200	6.5 mm/16	Washer for Dia. 6.5 mm Cancellous Screws	6
				4-184		Forceps	1

Basic Instrument Set for Screws and Plates in Graphic Case

Listing: T-340 **Basic Instrument Set for Screws and Plates**
(Steel Plates in Mat. DIN ISO 5832-1)



T-341 Tray for Basic Instrument Set (Steel)



19E-100.01 Standard Container, silver, 580 x 280 x 100 mm
(recommended)

Basic Instrument Set for Screws and Plates in Graphic Case

Listing:	T-340	Basic Instrument Set for Screws and Plates (Steel Plates in Mat. DIN ISO 5832-1)
	T-341	Tray for Basic Instrument Set (Steel)

Cat. No.:		Pcs.:
9-020	Drill Bit , \varnothing 3.2 mm , with end to fit quick coupling	3
9-026	Drill Bit , \varnothing 4.5 mm , with end to fit quick coupling	2
9-104	Malleolar Countersink , 3.2 \varnothing mm tip , for Malleolar Screws	1
9-092	Tap Handle for \varnothing 4.5 mm and \varnothing 6.5 mm Taps and Countersink	1
9-078	Tap , \varnothing 4.5 mm	2
9-079	Tap , \varnothing 6.5 mm	1
9-256	Tap Sleeve , \varnothing 3.5 mm (also for use as \varnothing 3.2 mm Drill Guide)	1
9-280	Insert Drill Sleeve , \varnothing 4.5 mm / \varnothing 3.2 mm , 58 mm Length	1
9-260	Tap Sleeve , \varnothing 4.5 mm (also for use as \varnothing 4.5 mm Drill Guide)	1
9-144	Hexagonal Screwdriver , width across flats 3.5 mm	1
9-176	Hexagonal Screwdriver for quick coupling , width across flats 3.5 mm	1
9-114	Depth Gauge for \varnothing 4.5 mm and \varnothing 6.5 mm Screws	1
9-116	Sharp Hook	1
9-262	Drill Sleeve for Plates , 40 mm Length , for Round Hole Plates	1
9-440	Articulated Tension Device with Gauge , 20 mm span	1
9-432	Socket Wrench , width across flats 11 mm	1
9-438	Combination Wrench , width across flats 11 mm	1
9-272	Neutral and Load Drill Guide , \varnothing 4.5 mm	1
3-070	Bending Template , 5 Holes	1
3-071	Bending Template , 7 Holes	1
3-072	Bending Template , 9 Holes	1
9-106	Countersink \varnothing 4.5 mm	1
7-264	Kirschner Wires \varnothing 2.0 mm / 150 mm (Steel)	1

Small Fragment Instrument and Steel Implant Set in Graphic Case

Listing: T-415 **Small Fragment Instrument and Steel Implant Set**
(Steel Screws and Plates in Mat. DIN ISO 5832-1)



T-416 Tray for Small Fragment Instrument and Steel Implant Set



19E-101.01 1/1 Standard Container, silver, 580 x 280 x 135 mm
(recommended)

Small Fragment Instrument and Steel Implant Set in Graphic Case

- Listing: T-415 Small Fragment Instrument and Steel Implant Set
(Steel Screws and Plates in Mat. DIN ISO 5832-1)**
- T-416 Tray for Small Fragment Instrument and Steel Implant Set

Cat.No.:	Instruments for Dia. 3.5mm / Dia. 4.0mm Screws	Pcs.:	Screw type used in (Dia. mm)
9-016	Drill Bit Dia. 2.5 mm with quick coupling end	2	4.0 mm (Cancellous Screws)
9-022	Drill Bit Dia. 3.5 mm with quick coupling end	2	3.5 mm (Cortex Screws)
9-102	Small Countersink with quick coupling end	1	
9-074	Cortical Tap Dia. 3.5 mm with quick coupling end	2	3.5 mm (Cortex Screws)
9-076	Cancellous Tap Dia. 4.0 mm with quick coupling end	2	4.0 mm (Cancellous Screws)
9-094	Handle with quick coupling	1	
9-252	Drill Guide and Drill Sleeve for Drill Bit Dia. 2.0 mm	1	
9-254	Tap Sleeve Dia.3.5 mm, Sleeve for Drill Bit Dia.2.5mm	1	
9-278	Insert Sleeve Dia.3.5mm / Dia. 2.5mm	1	
9-172	Small Hexagonal Screwdriver, insert for Hex.-Socket 2.5 mm	1	
9-174	Small Hexagonal Screwdriver Shaft, Hex. 2.5mm	1	
9-112	Small Depth Gauge	1	
9-270	Neutral Drill Guide Dia. 3.5 mm, used for Drill Bit 2.5 mm	1	

Cat.No.:	General Instruments	Pcs.:	Screw type used in (Dia. mm)
9-116	Sharp Hook	1	
9-402	Bending Iron for Small Plates	1	
9-404	Bending Iron for Small Plates	1	
3-075	Bending Template for Small Plates 7 Holes	1	
3-076	Bending Template for Small Plates 9 Holes	1	
9-580	Self-centring Bone Holding Forceps	1	
9-596	Reduction Forceps with points	1	
9-600	Reduction Forceps for Small Fragments	1	
9-650	Small Hohmann Retractor 6 mm	1	
9-652	Small Hohmann Retractor 8 mm	1	
9-654	Hohmann Retractor with broad shank	1	
9-758	Periosteal Elevator, width 6 mm small, straight edge	1	
4-184	Screw Forceps self holding	1	

Cat.No.:	Dia. 3.5 mm Cortex Screws	Pcs.:	Cat.No. :	Dia. 4.0 mm Cancellous Screws short Thread	Pcs.:
4-002	Length 10 mm , Steel	8	4-102	Length 12 mm / 5 mm, Steel	6
4-004	Length 12 mm , Steel	8	4-104	Length 14 mm / 5 mm, Steel	6
4-006	Length 14 mm , Steel	8	4-106	Length 16 mm / 6 mm, Steel	6
4-008	Length 16 mm , Steel	8	4-108	Length 18 mm / 7 mm, Steel	6
4-010	Length 18 mm , Steel	8	4-110	Length 20 mm / 8 mm, Steel	6
4-012	Length 20 mm , Steel	8	4-112	Length 22 mm / 9 mm, Steel	6
4-014	Length 22 mm , Steel	8	4-114	Length 24 mm / 10 mm, Steel	6
4-016	Length 24 mm , Steel	8	4-116	Length 26 mm / 12 mm, Steel	6
4-018	Length 26 mm , Steel	8	4-118	Length 28 mm / 14 mm, Steel	6
4-020	Length 28 mm , Steel	8	4-120	Length 30 mm / 14 mm, Steel	6
4-024	Length 32 mm , Steel	8	4-122	Length 35 mm / 14 mm, Steel	6
4-026	Length 36 mm , Steel	8	4-124	Length 40 mm / 14 mm, Steel	6
4-028	Length 40 mm , Steel	8	4-126	Length 45 mm / 15 mm, Steel	6
4-030	Length 45 mm , Steel	8	4-128	Length 50 mm / 15 mm, Steel	6
4-032	Length 50 mm , Steel	8			
			4-176	Small Washer for Dia. 3.5 mm / 4.0 mm Screws	6
			4-240	Small T-Plate right angle 3/3 Holes Steel	1
			4-244	Small T-Plate right angle 4/4 Holes Steel	1
			4-250	Small T-Plate oblique angle 3/3 Holes Steel	1
			4-252	Small T-Plate oblique angle 5/3 Holes Steel	1
			4-220	One Third Tubular Plate 2 Holes, Steel	1
			4-222	One Third Tubular Plate 3 Holes, Steel	1
			4-224	One Third Tubular Plate 4 Holes, Steel	1
			4-226	One Third Tubular Plate 5 Holes, Steel	1
			4-228	One Third Tubular Plate 6 Holes, Steel	1
			4-230	One Third Tubular Plate 7 Holes, Steel	1
			4-232	One Third Tubular Plate 8 Holes, Steel	1
			4-256	Cloverleaf Plate, Shaft 3 Holes, Steel	1
			4-258	Cloverleaf Plate, Shaft 4 Holes, Steel	1
			4-208	Small Fragment Plate 6 Holes, Steel	2
			4-210	Small Fragment Plate 7 Holes, Steel	2
			4-212	Small Fragment Plate 8 Holes, Steel	2
			4-214	Small Fragment Plate 10 Holes, Steel	2
			7-254	Kirschner Wire Dia.1.2 mm x Length 150 mm, Steel	12
			7-260	Kirschner Wire Dia.1.6 mm x Length 150 mm, Steel	12
			7-264	Kirschner Wire Dia.2.0 mm x Length 150 mm, Steel	12

Small Fragment Instrument and Titanium Implant Set in Graphic Case

Listing: T-4150 **Small Fragment Instrument and Titanium Implant Set
(Titanium Screws and Plates in Mat. DIN ISO 5832-3)**



T-416 Tray for Small Fragment Instrument and Titanium Implant Set



19E-101.01 1/1 Standard Container, silver, 580 x 280 x 135 mm
(recommended)

Small Fragment Instrument and Titanium Implant Set in Graphic Case

**Listing: T-4150 Small Fragment Instrument and Titanium Implant Set
(Titanium Screws and Plates in Mat. DIN ISO 5832-3)**

T-416 Tray for Small Fragment Instrument and Titanium Implant Set


Cat.No.:	Instruments for Dia. 3.5mm / Dia. 4.0mm Screws	Pcs.:	Screw type used in (Dia. mm)
9-016	Drill Bit Dia. 2.5 mm with quick coupling end	2	4.0 mm (Cancellous Screws)
9-022	Drill Bit Dia. 3.5 mm with quick coupling end	2	3.5 mm (Cortex Screws)
9-102	Small Countersink with quick coupling end	1	
9-074	Cortical Tap Dia. 3.5 mm with quick coupling end	2	3.5 mm (Cortex Screws)
9-076	Cancellous Tap Dia. 4.0 mm with quick coupling end	2	4.0 mm (Cancellous Screws)
9-094	Handle with quick coupling	1	
9-252	Drill Guide and Drill Sleeve for Drill Bit Dia. 2.0 mm	1	
9-254	Tap Sleeve Dia.3.5 mm, Sleeve for Drill Bit Dia.2.5mm	1	
9-278	Insert Sleeve Dia.3.5mm / Dia. 2.5mm	1	
9-172	Small Hexagonal Screwdriver, insert for Hex.-Socket 2.5 mm	1	
9-174	Small Hexagonal Screwdriver Shaft, Hex. 2.5mm	1	
9-112	Small Depth Gauge	1	
9-270	Neutral Drill Guide Dia. 3.5 mm, used for Drill Bit 2.5 mm	1	

Cat.No.:	General Instruments	Pcs.:	Screw type used in (Dia. mm)
9-116	Sharp Hook	1	
9-402	Bending Iron for Small Plates	1	
9-404	Bending Iron for Small Plates	1	
3-075	Bending Template for Small Plates 7 Holes	1	
3-076	Bending Template for Small Plates 9 Holes	1	
9-580	Self-centring Bone Holding Forceps	1	
9-596	Reduction Forceps with points	1	
9-600	Reduction Forceps for Small Fragments	1	
9-650	Small Hohmann Retractor 6 mm	1	
9-652	Small Hohmann Retractor 8 mm	1	
9-654	Hohmann Retractor with broad shank	1	
9-758	Periosteal Elevator, width 6 mm small, straight edge	1	
4-184	Screw Forceps self holding	1	

Cat.No.:	Dia. 3.5 mm Cortex Screws	Pcs.:	Cat.No. :	Dia. 4.0 mm Cancellous Screws short Thread	Pcs.:
4-00200	Length 10 mm , Titanium	8	4-10200	Length 12 mm / 5 mm, Titanium	6
4-00400	Length 12 mm , Titanium	8	4-10400	Length 14 mm / 5 mm, Titanium	6
4-00600	Length 14 mm , Titanium	8	4-10600	Length 16 mm / 6 mm, Titanium	6
4-00800	Length 16 mm , Titanium	8	4-10800	Length 18 mm / 7 mm, Titanium	6
4-01000	Length 18 mm , Titanium	8	4-11000	Length 20 mm / 8 mm, Titanium	6
4-01200	Length 20 mm , Titanium	8	4-11200	Length 22 mm / 9 mm, Titanium	6
4-01400	Length 22 mm , Titanium	8	4-11400	Length 24 mm / 10 mm, Titanium	6
4-01600	Length 24 mm , Titanium	8	4-11600	Length 26 mm / 12 mm, Titanium	6
4-01800	Length 26 mm , Titanium	8	4-11800	Length 28 mm / 14 mm, Titanium	6
4-02000	Length 28 mm , Titanium	8	4-12000	Length 30 mm / 14 mm, Titanium	6
4-02400	Length 32 mm , Titanium	8	4-12200	Length 35 mm / 14 mm, Titanium	6
4-02600	Length 36 mm , Titanium	8	4-12400	Length 40 mm / 14 mm, Titanium	6
4-02800	Length 40 mm , Titanium	8	4-12600	Length 45 mm / 15 mm, Titanium	6
4-03000	Length 45 mm , Titanium	8	4-12800	Length 50 mm / 15 mm, Titanium	6
4-03200	Length 50 mm , Titanium	8			
			4-1760	Small Washer for Dia. 3.5 mm / 4.0 mm Screws	6
			4-2400	Small T-Plate right angle 3/3 Holes Titanium	1
			4-2440	Small T-Plate right angle 4/4 Holes Titanium	1
			4-2500	Small T-Plate oblique angle 3/3 Holes Titanium	1
			4-2520	Small T-Plate oblique angle 5/3 Holes Titanium	1
			4-2200	One Third Tubular Plate 2 Holes, Titanium	1
			4-2220	One Third Tubular Plate 3 Holes, Titanium	1
			4-2240	One Third Tubular Plate 4 Holes, Titanium	1
			4-2260	One Third Tubular Plate 5 Holes, Titanium	1
			4-2280	One Third Tubular Plate 6 Holes, Titanium	1
			4-2300	One Third Tubular Plate 7 Holes, Titanium	1
			4-2320	One Third Tubular Plate 8 Holes, Titanium	1
			4-2560	Cloverleaf Plate, Shaft 3 Holes, Titanium	1
			4-2580	Cloverleaf Plate, Shaft 4 Holes, Titanium	1
			4-2080	Small Fragment Plate 6 Holes, Titanium	2
			4-2100	Small Fragment Plate 7 Holes, Titanium	2
			4-2120	Small Fragment Plate 8 Holes, Titanium	2
			4-2140	Small Fragment Plate 10 Holes, Titanium	2
			7-254	Kirschner Wire Dia.1.2 mm x Length 150 mm	12
			7-260	Kirschner Wire Dia.1.6 mm x Length 150 mm	12
			7-264	Kirschner Wire Dia.2.0 mm x Length 150 mm	12

Mini Instrument and Steel Implant Set in Graphic Case

**Listing: T-420 Mini Instrument and Steel Implant Set
(Steel Screws and Plates in Mat. DIN ISO 5832-1)**



T-421 Tray for Mini Instrument and Implant Set (Steel)



19E-100.01 Standard Container, silver, 580 x 280 x 100 mm
(recommended)

Mini Instrument and Steel Implant Set in Graphic Case

**Listing: T-420 Mini Instrument and Steel Implant Set
(Steel Screws and Plates in Mat. DIN ISO 5832-1)**

T-421 Tray for Mini Instrument and Implant Set (Steel)

Cat.No.:	Instruments for Dia. 2.0mm / Dia. 2.7mm Screws	Pcs.:	Screw type used in (Dia. mm)
9-014	Drill Bit Dia. 2.0 mm with quick coupling	2	2.0 mm + 2.7 mm
9-018	Drill Bit Dia. 2.7 mm with quick coupling	2	2.7 mm
9-102	Small Countersink with end to fit quick coupling Handle	1	2.7 mm
9-073	Tap Dia. 2.7 mm with end to fit quick coupling Handle	2	2.7 mm
9-094	Handle with quick coupling	1	
9-252	Drill Guide and Drill Sleeve for Dia. 2.0 mm Drill Bit	1	2.0 mm + 2.7 mm
9-256	Tap Sleeve Dia. 3.5 mm	1	2.7 mm
9-172	Small Hexagonal Screwdriver 2.5 mm with holding Sleeve	1	2.7 mm
9-112	Depth Gauge (only for 2.7 mm Screws)	1	2.7 mm

Cat.No.:	Instruments for Dia. 1.5 mm / Dia. 2.0 mm Mini Screws	Pcs.:	Screw type used in (Dia. mm)
9-010	Drill Bit Dia. 1.1 mm with quick coupling	2	1.5 mm
9-012	Drill Bit Dia. 1.5 mm with quick coupling	2	2.0 mm
9-100	Mini Countersink to fit Handle with mini quick coupling	1	1.5 mm + 2.0 mm
9-071	Tap Dia. 1.5 mm to fit Handle with mini quick coupling	2	1.5 mm
9-072	Tap Dia. 2.0 mm to fit Handle with mini quick coupling	2	2.0 mm
9-096	Handle with mini quick coupling	1	
9-250	Mini Drill Sleeve for Dia. 1.1 and Dia. 1.5 mm Drill Bits	1	1.5 mm + 2.0 mm
9-110	Mini Depth Gauge (only for 2.0 mm Screws)	1	2.0 mm

Cat.No.:	General Instruments	Pcs.:
9-116	Sharp Hook	1
9-170	Hexagonal Screw Driver 1.5 mm with holding Sleeve	1
9-400	Bending Iron for Finger Plates1	
9-630	Holding Forceps for small Plates	1
9-596	Reduction Forceps with points	1
9-650	Small Retractor, 6 mm	1
9-652	Small Retractor, 8 mm	1
9-654	Retractor for small fragments with broad shank , short	2
9-760	Periosteal Elevator, small with straight edge 3 mm	1

Cat.No.:	Dia. 2.7 mm Cortex Screws	Pcs.:	Cat.No. :	Dia. 1.5 mm Cortex Screws	Pcs.:
4-460	Length 6 mm , Steel	6	4-400	Length 6 mm , Steel	4
4-462	Length 8 mm , Steel	6	4-402	Length 7 mm , Steel	4
4-464	Length 10 mm , Steel	6	4-404	Length 8 mm , Steel	4
4-466	Length 12 mm , Steel	6	4-406	Length 9 mm , Steel	4
4-468	Length 14 mm , Steel	6	4-408	Length 10 mm , Steel	4
4-470	Length 16 mm , Steel	6	4-410	Length 11 mm , Steel	4
4-472	Length 18 mm , Steel	6	4-412	Length 12 mm , Steel	4
4-474	Length 20 mm , Steel	6	4-414	Length 14 mm , Steel	4
4-476	Length 22 mm , Steel	6	4-416	Length 16 mm , Steel	4
4-478	Length 24 mm , Steel	6			


Cat.No.:	Dia. 2.0 mm Cortex Screws	Pcs.:
4-430	Length 6 mm , Steel	6
4-432	Length 8 mm , Steel	6
4-434	Length 10 mm , Steel	6
4-436	Length 12 mm , Steel	6
4-438	Length 14 mm , Steel	6
4-440	Length 16 mm , Steel	6
4-442	Length 18 mm , Steel	6
4-444	Length 20 mm , Steel	6

Cat.No.:	Plates for Dia. 2.7 mm Cortex Screws	Pcs.:	Cat.No.:	Plates for Dia. 1.5 mm / 2.0 mm Cortex Screws	Pcs.:
4-540	Quarter Tubular Plate Length 25 mm / 3 Hole , Steel	2	4-590	Mini Plate 2.0 straight , 3 Hole / Length 17 mm , Steel	2
4-542	Quarter Tubular Plate Length 33 mm / 4 Hole , Steel	2	4-592	Mini Plate 2.0 straight , 4 Hole / Length 23 mm , Steel	2
4-544	Quarter Tubular Plate Length 41 mm / 5 Hole , Steel	2	4-594	Mini Plate 2.0 straight , 5 Hole / Length 29 mm , Steel	2
4-546	Quarter Tubular Plate Length 49 mm / 6 Hole , Steel	2	4-600	Mini L - Plate 2.0 oblique left 2/2 Hole , Steel	2
4-548	Quarter Tubular Plate Length 57 mm / 7 Hole , Steel	2	4-602	Mini L - Plate 2.0 oblique right 2/2 Hole , Steel	2
4-550	Quarter Tubular Plate Length 65 mm / 8 Hole , Steel	2	4-604	Mini T - Plate 2.0 , 2/2 Hole , Steel	4
4-552	L - Plate , oblique left 3/2 Hole , Steel	2			
4-554	L - Plate , oblique right 3/2 Hole , Steel	2			
4-556	T - Plate , Shaft 3 Hole / Head 2 Hole , Steel	4			

Cat.No.:	Kirschner Wires with Trocar	Pcs.:
7-200	Dia. 0.8 mm / Length 70 mm , Steel	10
7-201	Dia. 1.0 mm / Length 70 mm , Steel	10
7-252	Dia. 1.0 mm / Length 150 mm , Steel	10
7-254	Dia. 1.2 mm / Length 150 mm , Steel	10
7-260	Dia. 1.6 mm / Length 150 mm , Steel	10

Mini Instrument / Titanium Implant Set in Graphic Case

Listing: T-4200 Mini Instrument / Titanium Implant Set
(Titanium Screws and Plates in Mat. DIN ISO 5832-3)



T-421 Tray for Mini Instrument and Implant Set (Steel)



19E-100.01 Standard Container, silver, 580 x 280 x 100 mm
(recommended)

Mini Instrument / Titanium Implant Set in Graphic Case

**Listing: T-4200 Mini Instrument / Titanium Implant Set
(Titanium Screws and Plates in Mat. DIN ISO 5832-3)**

T-421 Tray for Mini Instrument and Implant Set (Steel)

Cat.No.:	Instruments for Dia. 2.0mm / Dia. 2.7mm Screws	Pcs.:	Screw type used in (Dia. mm)
9-014	Drill Bit Dia. 2.0 mm with quick coupling	2	2.0 mm + 2.7 mm
9-018	Drill Bit Dia. 2.7 mm with quick coupling	2	2.7 mm
9-102	Small Countersink with end to fit quick coupling Handle	1	2.7 mm
9-073	Tap Dia. 2.7 mm with end to fit quick coupling Handle	2	2.7 mm
9-094	Handle with quick coupling	1	
9-252	Drill Guide and Drill Sleeve for Dia. 2.0 mm Drill Bit	1	2.0 mm + 2.7 mm
9-256	Tap Sleeve Dia. 3.5 mm	1	2.7 mm
9-172	Small Hexagonal Screwdriver 2.5 mm with holding Sleeve	1	2.7 mm
9-112	Depth Gauge (only for 2.7 mm Screws)	1	2.7 mm

Cat.No.:	Instruments for Dia. 1.5 mm / Dia. 2.0 mm Mini Screws	Pcs.:	Screw type used in (Dia. mm)
9-010	Drill Bit Dia. 1.1 mm with quick coupling	2	1.5 mm
9-012	Drill Bit Dia. 1.5 mm with quick coupling	2	2.0 mm
9-100	Mini Countersink to fit Handle with mini quick coupling	1	1.5 mm + 2.0 mm
9-071	Tap Dia. 1.5 mm to fit Handle with mini quick coupling	2	1.5 mm
9-072	Tap Dia. 2.0 mm to fit Handle with mini quick coupling	2	2.0 mm
9-096	Handle with mini quick coupling	1	
9-250	Mini Drill Sleeve for Dia. 1.1 and Dia. 1.5 mm Drill Bits	1	1.5 mm + 2.0 mm
9-110	Mini Depth Gauge (only for 2.0 mm Screws)	1	2.0 mm

Cat.No.:	General Instruments	Pcs.:
9-116	Sharp Hook	1
9-170	Hexagonal Screw Driver 1.5 mm with holding Sleeve	1
9-400	Bending Iron for Finger Plates1	
9-630	Holding Forceps for small Plates	1
9-596	Reduction Forceps with points	1
9-650	Small Retractor, 6 mm	1
9-652	Small Retractor, 8 mm	1
9-654	Retractor for small fragments with broad shank , short	2
9-760	Periosteal Elevator, small with straight edge 3 mm	1

Cat.No.:	Dia. 2.7 mm Cortex Screws	Pcs.:	Cat.No. :	Dia. 1.5 mm Cortex Screws	Pcs.:
4-46000	Length 6 mm , Ti	6	4-40000	Length 6 mm , Ti	4
4-46200	Length 8 mm , Ti	6	4-40200	Length 7 mm , Ti	4
4-46400	Length 10 mm , Ti	6	4-40400	Length 8 mm , Ti	4
4-46600	Length 12 mm , Ti	6	4-40600	Length 9 mm , Ti	4
4-46800	Length 14 mm , Ti	6	4-40800	Length 10 mm , Ti	4
4-47000	Length 16 mm , Ti	6	4-41000	Length 11 mm , Ti	4
4-47200	Length 18 mm , Ti	6	4-41200	Length 12 mm , Ti	4
4-47400	Length 20 mm , Ti	6	4-41400	Length 14 mm , Ti	4
4-47600	Length 22 mm , Ti	6	4-41600	Length 16 mm , Ti	4
4-47800	Length 24 mm , Ti	6			

Cat.No.:	Dia. 2.0 mm Cortex Screws	Pcs.:
4-43000	Length 6 mm , Ti	6
4-43200	Length 8 mm , Ti	6
4-43400	Length 10 mm , Ti	6
4-43600	Length 12 mm , Ti	6
4-43800	Length 14 mm , Ti	6
4-44000	Length 16 mm , Ti	6
4-44200	Length 18 mm , Ti	6
4-44400	Length 20 mm , Ti	6



Cat.No.:	Plates for Dia. 2.7 mm Cortex Screws	Pcs.:
4-5400	Quarter Tubular Plate Length 25 mm / 3 Hole , Ti	2
4-5420	Quarter Tubular Plate Length 33 mm / 4 Hole , Ti	2
4-5440	Quarter Tubular Plate Length 41 mm / 5 Hole , Ti	2
4-5460	Quarter Tubular Plate Length 49 mm / 6 Hole , Ti	2
4-5480	Quarter Tubular Plate Length 57 mm / 7 Hole , Ti	2
4-5500	Quarter Tubular Plate Length 65 mm / 8 Hole , Ti	2
4-5520	L - Plate , oblique left 3/2 Hole , Ti	2
4-5540	L - Plate , oblique right 3/2 Hole , Ti	2
4-5560	T - Plate , Shaft 3 Hole / Head 2 Hole , Ti	4

Cat.No.:	Plates for Dia. 1.5 mm / 2.0 mm Cortex Screws	Pcs.:
4-5900	Mini Plate 2.0 straight , 3 Hole / Length 17 mm , Ti	2
4-5920	Mini Plate 2.0 straight , 4 Hole / Length 23 mm , Ti	2
4-5940	Mini Plate 2.0 straight , 5 Hole / Length 29 mm , Ti	2
4-6000	Mini L - Plate 2.0 oblique left 2/2 Hole , Ti	2
4-6020	Mini L - Plate 2.0 oblique right 2/2 Hole , Ti	2
4-6040	Mini T - Plate 2.0 , 2/2 Hole , Ti	4

Cat.No.:	Kirschner Wires with Trocar	Pcs.:
7-200	Dia. 0.8 mm / Length 70 mm , Steel	10
7-201	Dia. 1.0 mm / Length 70 mm , Steel	10
7-252	Dia. 1.0 mm / Length 150 mm , Steel	10
7-254	Dia. 1.2 mm / Length 150 mm , Steel	10
7-260	Dia. 1.6 mm / Length 150 mm , Steel	10

Implant Set for Intramedullary Gliding Nails

Listing: **T-5000** **Implant Set for Intramedullary Gliding Nails**
(Steel and Titanium Gliding Nail Set)



1-5001 Tray (empty)



19E-100.01 Standard Container, silver, 580 x 280 x 100 mm
(recommended)

Implant Set for Intramedullary Gliding Nails

**Listing: T-5000 Implant Set for Intramedullary Gliding Nails
(Steel and Titanium Gliding Nail Set)**

1-5001 Tray (empty)

Steel Intramedullary Gliding Nail Implants (Stainless Steel DIN ISO 5832-1)


1-5250	Dia. 1,5 mm x Length 150 mm	4 Pcs.
1-5251	Dia. 1,5 mm x Length 400 mm	4 Pcs.
1-5252	Dia. 2,0 mm x Length 200 mm	4 Pcs.
1-5253	Dia. 2,0 mm x Length 400 mm	4 Pcs.
1-5254	Dia. 2,5 mm x Length 250 mm	4 Pcs.
1-5255	Dia. 2,5 mm x Length 400 mm	4 Pcs.
1-5256	Dia. 3,0 mm x Length 300 mm	4 Pcs.
1-5257	Dia. 3,0 mm x Length 400 mm	4 Pcs.
1-5258	Dia. 3,5 mm x Length 350 mm	4 Pcs.
1-5259	Dia. 3,5 mm x Length 400 mm	4 Pcs.

Titanium Intramedullary Gliding Nail Implants (Titanium DIN ISO 5832-3)

1-5263	Dia. 2,0 mm x Length 450 mm Colour Coding Green	4 Pcs.
1-5264	Dia. 2,5 mm x Length 450 mm Colour Coding Pink	4 Pcs.
1-5265	Dia. 3,0 mm x Length 450 mm Colour Coding Gold	4 Pcs.
1-5266	Dia. 3,5 mm x Length 450 mm Colour Coding Blue	4 Pcs.
1-5267	Dia. 4,0 mm x Length 450 mm Colour Coding Violet	4 Pcs.
1-5268	Dia. 4,5 mm x Length 450 mm Colour Coding Gray	4 Pcs.
1-5269	Dia. 5,0 mm x Length 450 mm Colour Coding Silver	4 Pcs.

Standard Set External Fixateur (Tubular System) in Graphic Case

Listing: T-8000 Grafic Case Standard Set External Fixateur (Tubular System)



1-8001 Tray (empty)



19E-100.01 Standard Container, silver, 580 x 280 x 100 mm
(recommended)

Standard Set External Fixateur (Tubular System) in Graphic Case

Listing: T-8000 Grafic Case Standard Set External Fixateur (Tubular System)

1-8001 Tray (empty)

Cat.No.		Pcs. :
9-023	Drill Bit Ø 3.5 mm , extra long	2
9-438	Combination Wrench with across flats 11 mm	1
8-0068	Simple T-Handle for Steinmann Pins and Schanz Screws Ø5.0 mm	1
8-0070	Socket Wrench with across flats 11mm 180mm long for threadet conical Bolts and External Fixateur	1
8-0072	Hexagonal Allen Key withacross flats 3.5 mm	1
8-0092	Trocar with point Ø 3.5 mm short	1
8-0061	Drill Sleeve Ø 5.0 mm / Ø 3.5 mm short	1
8-0065	Drill Sleeve Ø 6.0 mm / Ø 5.0 mm short	1
8-0062	Trocar with point for Aiming Device Ø 3.5 mm long	1
8-0060	Drill Sleeve Ø 5.0 mm / Ø 3.5 mm long	1
8-0064	Drill Sleeve Ø 6.0 mm / Ø 5.0 mm long	1

Cat.No.	Fixateur Components	Pcs. :
8-0080	Open Compressor	2
8-0100	Tube Rod Ø 11 mm , 100 mm Length (Steel)	1
8-0102	Tube Rod Ø 11 mm , 125 mm Length (Steel)	1
8-0104	Tube Rod Ø 11 mm , 150 mm Length (Steel)	2
8-0106	Tube Rod Ø 11 mm , 200 mm Length (Steel)	2
8-0108	Tube Rod Ø 11 mm , 250 mm Length (Steel)	2
8-0110	Tube Rod Ø 11 mm , 300 mm Length (Steel)	2
8-0112	Tube Rod Ø 11 mm , 350 mm Length (Steel)	2
8-0114	Tube Rod Ø 11 mm , 400 mm Length (Steel)	2
8-0170	Plug for Tube Rods	10
8-0188	Open adjustable Clamp	12
8-0190	Ø 11 mm Tube-to Tube Clamp	2

**Schanz Screws Ø 5.0 mm , Thread - Length 50 mm (with Trocar Point and Triangular Shaft)
(Steel Implants DIN ISO 5832-1)**

Cat.No.		Pcs. :
8-1264	Ø 5.0 mm , Length 125 mm	6
8-1266	Ø 5.0 mm , Length 150 mm	6
8-1268	Ø 5.0 mm , Length 170 mm	6
8-1270	Ø 5.0 mm , Length 200 mm	6



Technische Änderungen vorbehalten
Technical features are subject to change
Cambios técnicos reservados
Ci riserviamo il diritto di modifiche tecniche



erbrich - instrumente GmbH
Eugenstr. 33 · D-78532 Tuttlingen
Tel.: +49 74 61 / 96 57 30
Fax: +49 74 61 / 9657 320
E-mail: info@erbrich-instrumente.de
Website: www.erbrich-instruments.com